



Vireo Monitor Assembly Instructions

May 2025



Unpacking

Open the box from the top

1. Cables
2. Protective foam

To access monitor, remove cables and protective foam.



3 mechanical components



1. Monitor base
2. Monitor stand
3. 32 inch Monitor





Connecting the monitor

When installing the arm, insert the upper portion first, then press on use the latch button for smooth assembly.





Cable Connectivity

4 Cables to connect:

1. Power
2. HDMI
3. USB
4. Keyboard and Mouse Dongle





Disassembly and packing

1. Remove all the cables.
2. Disconnect the arm by pressing on the latch button to release the monitor stand.

Place components as shown in the unpacking instructions

Questions:

Email us at help@ramonaoptics.com

